

Laser Hair Plus User Manual

1. The Laser Hair Plus has 8 independent Laser beams, disperse laser with patented scattering devices.
2. Insert the power adapter attachment into the handle end of the Laser Hair Plus. Plug the power adaptor into any convenient wall electrical outlet.
3. For the initial charge, the LED indicator will blink red. Strongly recommend charge the battery at least 4 hours, and the LED indicator will glow red when fully charged.
4. After fully charge, the LED indicator (red color) will go out (not be lit) when power adapter be reconnected.
5. For safety reason, the laser device will be cut-off automatically while battery is charging, without any function to switch the ON/OFF button.(Can hear one beep tone and blue LED blink)
6. Plug out the charging cable, and push the ON/OFF button to turn on the Laser Hair Plus. The LED indicator will light blue. You can clearly see the laser light when the device is operating.
7. An audible beeping tone will be heard every 4 seconds with this activation of the Laser Hair Plus. The tone is added to help you determine the time period between movements. Hold the comb on one spot for 4 seconds, then move it approximately 1/3 inch (half distance between two rows of comb teeth) to the next spot
8. Move the Laser Hair Plus against the direction of the hair growth – front to back and bottom to top.
9. The Laser Hair Plus intelligently manages treatment time in 20 minutes session. Use a powerful onboard microprocessor system that is capable of managing exposure session in each treatment, the laser device will be turn-off automatically after 20 minutes. Needless to say, this is a complicated procedure for most home phototherapy patients and accidental over exposure would often result.
10. Recommended usage of the Laser Hair Plus is one programming time (20 minutes) for risky hair loss patients, and two programming time (40 minutes) for bald hair loss patients, three treatments a week. There is no evidence to suggest over exposure is possible or side effects exist. However, it is recommended to abide by the suggested exposure times. Longer exposure times may not produce better results from treatment.
11. The LED indicator should blink blue in 60 seconds when the unit indicates low battery.
12. If the low battery signal is left unattended in a period, Laser Hair Plus may return to SLEEP mode, and you can not find any response when push the ON/OFF button. Perform the WAKE function for any time you find the situation, to recharge the battery again will active the WAKE mode.
13. Follow step 2 to plug the power adapter cable without push the ON/OFF button. You will find all lasers/ red LED/ blue LED/ beep alert in first second, which means the WAKE mode is turning on.



CAUTION: LASER LIGHT WHEN OPEN. DO NOT STARE INTO BEAM. KEEP OUT OF REACH OF ALL CHILDREN